

**LIQUIDS**



**PASTES**

**INSULATION**

# **THW-L2**


Conforms to standard ASTM D7896-14

*Portable Thermal Conductivity  
Meter for Measurement of  
Liquids, Pastes, & Insulation.*

# **THERMTEST PORTABLE**



# PORTABLE. ECONOMICAL. ACCURATE.




## STANDARD THW-L2 SENSOR


The THW-L2 comes with the standard 60 mm sensor for testing of liquids and pastes. The sensor is fully inserted into an isothermal sample and measurements are made with the push of a button. In less than 2 seconds, results are displayed for thermal conductivity. Volumes as small as 15 mL can be tested. The easy to use THW software enables measurements to be controlled via a computer. When used in conjunction with the optional thermoelectric dry bath, the software can even automate measurements of thermal conductivity with temperature.

## THERMAL CONDUCTIVITY OF INSULATION

The THW-L2 is capable of testing low thermal conductivity, compressible insulations, without the liquid sample holder. Samples with minimum dimensions of 65 mm in length and 2 mm in thickness can be tested in a horizontal configuration, with the sensor sandwiched between the samples. Using the THW-L2, expanded polystyrene was measured as  $0.031 \text{ W/m}\cdot\text{K}$ , which is within 5% of the literature values.



# SAMPLE MEASUREMENT




## THE SAMPLE

Using the THW-L2, liquid samples of unlimited size can be tested. With the small volume cell insert, volumes as small as 15 mL can be measured. The instrument's ability to control convection, using short test times, allows accurate measurements of a wide range of samples, with varying viscosities.



1 min.



## INSERT SENSOR

Once the liquid has been poured into the sample holder, the sensor can be vertically inserted in place. For smaller sample volumes, position the small volume cell into the sample holder, prior to the addition of the sample and sensor.



< 1 min.

# EFFICIENCY WITH EASE





TEMPSTEPWINDOW	
Sample Name	DIUF
Test Power (mW)	0
Test Time (s)	1.5
Start Temperature (c)	5
Temperature Step (c)	5
End Temperature (c)	70
Number of Tests	3
Time between tests (m)	5
<input type="button" value="run test"/> <input type="button" value="cancel"/>	

## RUN EXPERIMENT

The instrument can be controlled through the front panel with the push of a button, or through a Windows based computer software. Single temperature point measurements, as well as scheduled temperature intervals may be performed using the THW software.



< 1 min.



## EXPORTING RESULTS


Using the Windows based software, users can save, print, and export test results to Microsoft Excel, for further processing once testing is complete.



1 min.

## LIQUIDS VS. TEMPERATURE


Using the Thermoelectric Dry Bath, the Thermtest THW-L2 is capable of testing liquids over a temperature range of 0 to 100°C (32 to 212°F). In this application, deionized ultra-filtered (DIUF) water was tested from 10 to 70°C, using the THW-L2. All test results were within 5% of the literature values.



## OPTIONAL DRY BATH

The optional dry bath allows for automated measurements of thermal conductivity with temperature (-10 to 90°C or 0 to 100°C). The THW software automatically controls the increments of temperature steps between the desired start and end temperatures to ensure isothermal conditions and test results are free of temperature drift effects. For wider temperature ranges, additional cooling/heating options are available.





Thermtest has been advancing the measurement of thermal conductivity, thermal diffusivity, and specific heat for more than a decade. With more than 1500 satisfied customers, our unique combination of advanced thermal conductivity instrumentation for the laboratory, portable meters for the field, and accessories, enables us to provide ideal solutions to fit any material testing application and budget. Our proud commitment to being a leader in thermal conductivity has fueled our success through rigorous development and key partnerships, creating a lineup of industry leading testing solutions for the laboratory, field, and production-line.

